

Bioneer Corporation
8-11 Munpyeongseo-ro, Daedeok-gu, Daejeon, 34302
Republic of Korea
Tel: +82-42-930-8777 (Korea : 1588-9788)
Fax: +82-42-930-8688 E-mail: sales@bioneer.com

Bioneer Inc.
1301 Marina Village PKWY, Suite 110
Alameda, CA 94501, USA
Toll Free: +1-877-264-4300 Fax: +1-510-865-0350
E-mail: order.usa@bioneer.cus.com

Bioneer R&D Center
Korea Bio Park BLDG #B-702, 700 Daewangpangyo-ro
Bundang-gu, Seongnam-si, Gyeonggi-do, 13488
Republic of Korea
Tel: +82-31-628-0500 Fax: +82-31-628-0555

AccuPower[®] RT-PCR PreMix contains all the components necessary for cDNA synthesis and amplification in one tube.

Advantages

- Speed** Substantial reduction in reaction setup time.
- Stability** As each tube of *AccuPower[®]* RT-PCR PreMix contains a stabilizer (patented in US and Korea), which can maintain the stability of the RTase and DNA polymerase up to 2 years at -20°C. Pre-denaturation during PCR has no effects on the activity of the DNA polymerase.
- Reproducibility and Yield** The strict functional QC assays demonstrated highly consistent and reproducible RT/PCR performance. In most applications an increase in yield is observed as compared to the standard reactions.
- Simplicity** The fewer manual steps allow reduction in potential errors. Each tube contains tracking dye and precipitant for agarose gel electrophoresis, eliminating the needs for a separate loading buffer.

Experimental Protocol

- Mix the template RNA and the reverse primer in a sterile tube as indicated below;

Concentration of template RNA and primer			
Reaction volume		20 µl reaction	50 µl reaction
Template RNA	Total RNA	0.5~1.0 µg	1.0~2.0 µg
	Poly(A) RNA	0.05~0.1 µg	0.1~0.2 µg
Primer	Primer	10~30 pmole	20~50 pmole

- Incubate the mixture at 70°C for 5min and place it on ice.
- Transfer the incubated mixture and the forward primer to *AccuPower[®]* RT/PCR PreMix tube, then fill up the reaction volume with DEPC-DW.
- Dissolve the lyophilized blue pellet by vortexing, and briefly spin down.
- Add mineral oil to each tube (This step is unnecessary when using a thermal cycler with top heating.).
- Perform cDNA synthesis reaction as follows :
 - 42°C, 60 min. (cDNA synthesis)
 - 94°C, 5 min. (RTase inactivation)
- Perform PCR cycles according to the PCR condition.
(Annealing temperature and time need to be optimized for each primer/template combination.)

Trademark

AccuPower is a trademark of Bioneer Corporation.

Ordering Information

Tube type	Reaction	Cat.No	Description	Tube type	Reaction	Cat.No	Description
0.2ml Tube	20 ul	K-2055	0.2 ml thin-wall 8-strip tubes with attached cap / 96 tubes	96-well	20 ul	K-2262-4	thin-wall 96-well flat plate
		K-2055-B	0.2 ml thin-wall 8-strip tubes with attached cap / 480 tubes			K-2262-5	thin-wall 96-well full-skirted plate
	50 ul	K-2057	0.2 ml thin-wall 8-strip tubes with attached cap / 96 tubes			K-2262-6	thin-wall 96-well semi-skirted plate
		K-2057-B	0.2 ml thin-wall 8-strip tubes with attached cap / 480 tubes	384-well	5 ul	K-2084-1	thin-wall 384-well full-skirted plate
96-well	10 ul	K-2262-1	thin-wall 96-well flat plate		10 ul	K-2084-2	thin-wall 384-well full-skirted plate
		K-2262-2	thin-wall 96-well full-skirted plate		20 ul	K-2084-3	thin-wall 384-well full-skirted plate
		K-2262-3	thin-wall 96-well semi-skirted plate	0.5ml Tube	50ul	k-2056	0.5 ml thin-wall tubes with attached cap / 100 tubes

A complete product list appears on our web site at www.bioneer.com.