



## Zobell Marine Agar 2216

M384

### Intended use

Zobell Marine Agar 2216 is recommended for cultivation, isolation and enumeration of heterotrophic marine bacteria.

### Composition\*\*

Ingredients	Gms / Litre
Peptone	5.000
Yeast extract	1.000
Ferric citrate	0.100
Sodium chloride	19.450
Magnesium chloride	8.800
Sodium sulphate	3.240
Calcium chloride	1.800
Potassium chloride	0.550
Sodium bicarbonate	0.160
Potassium bromide	0.080
Strontium chloride	0.034
Boric acid	0.022
Sodium silicate	0.004
Sodium fluoride	0.0024
Ammonium nitrate	0.0016
Disodium phosphate	0.008
Agar	15.000
Final pH ( at 25°C)	7.6±0.2

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

Suspend 55.25 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and Pour in sterile Petri plates.

### Principle And Interpretation

Microorganisms in an aquatic environment may occur at all depths ranging from the surface region to the very bottom of the ocean trenches. The top layers and the bottom sediments harbour higher concentration of microorganisms (1). Marine microorganisms are vital to ecological cycles because they form the foundations of many food chains (2). Zobell Marine Agar formulated by Zobell (3), has a composition that mimics seawater (4) and thus helps the marine bacteria to grow abundantly. This medium has been used for the growth of marine bacteria (5, 6).

Zobell Marine Agar 2216 contains the nutrients, which are required for the growth of marine bacteria. These media have minerals as in seawater (7) and peptone and yeast extract as the sources of nutrients for the marine bacteria as reported by Jones (8). High amount of salt content is used to simulate seawater. Other minerals are used to mimic the mineral composition of seawater.

Pour plate and spread plate techniques can be used for enumeration. In the pour plate technique, the agar must be cooled to 42°C before inoculation to support thermo-sensitive nature of most marine bacteria. In spread plate technique, the medium is poured while still hot and allowed to cool and solidify before inoculation.

### Type of specimen

Marine water samples

### Specimen Collection and Handling:

For marine water samples follow appropriate techniques for handling specimens as per established guidelines ( ). After use, contaminated materials must be sterilized by autoclaving before discarding.

**Please refer disclaimer Overleaf.**

## Warning and Precautions :

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

## Limitations

1. The medium is recommended for the isolation of marine bacteria. Further biochemical and seological testing must be carried out for further identification.

## Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

## Quality Control

### Appearance

Cream to yellow homogeneous free flowing powder

### Gelling

Firm, comparable with 1.5% Agar gel.

### Colour and Clarity of prepared medium

Yellow coloured opalescent gel forms in Petri plates.

### Reaction

Reaction of 5.53% w/v aqueous solution at 25°C. pH : 7.6±0.2

### pH

7.40-7.80

### Cultural Response

M384: Cultural characteristics observed after an incubation at 20-25°C for 40-72 hours .

Organism	Inoculum (CFU)	Growth	Recovery
<b>Cultural Response</b>			
<i>Vibrio fischeri</i> ATCC 7744	50-100	good-luxuriant	≥50%
<i>Vibrio harveyi</i> ATCC 14126	50-100	good-luxuriant	≥50%

## Storage and Shelf Life

Store below 10-30°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period. .

## Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (9,10).

## Reference

1. Pelczar M.J., Jr., Reid R.D., Chan E.C.S., 1977, Microbiology, 4th Ed., Tata McGraw-Hill Publishing Company Ltd, New Delhi
2. Alcamo E.I., 2001, Fundamentals of Microbiology, 6th Ed., Jones AND Barlett Publishers
3. ZoBell C. E., 1941, J. Mar. Res., 4:42.
4. Lyman J. and Fleming R. H., 1940, J. Mar. Res. 3:134.
5. Sizemore R. K. and Stevenson L. H., 1970, Appl. Microbiol., 20:991.
6. Weiner R. M., Segall A. M. and Colwell R. R., 1985, Appl. Environ. Microbiol., 49:83.
7. Zobell C. E., 1940, J. Marine Research , 3:134
8. Jones, 1960, Bact. Proc. Pg. 36 (A29)

9. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2<sup>nd</sup> Edition.
10. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision : 02 / 2017

**Disclaimer :**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.